



Anti-Nef (HIV-1/CladeB)

Product Name	Anti- Nef (Clade B)
Catalog No	IT-001-014
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 Nef(B) (HXBc2) (Genebank No. K03455)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1 Nef(B). Cross-reactivity to Nef from other clades not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

