

# Anti-H1(California/06/2009)



Product Name	Anti-H1(California/06/2009(H1N1))
Catalog No	IT-003-001M4
Source	Mouse monoclonal antibody, clone 26-D11
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant HA1(A/California/06/2009(H1N1)) protein (a.a. 18-344)(Genebank accession <a href="#">ACP41935</a> ).
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H1(2009 H1N1, or swine flu), not react with seasonal strains (post 1957). No cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

