

Anti-NS4a (HCV)



Product Name	Anti- NS4a (HCV)
Catalog No	IT-004-018
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Peptide from a conserved region of NS4a polypeptide of Hepatitis C virus (a.a. 1055-1064) (RDKNQVEGEV).
Applications	WB(1:1000), ELISA(1:2000-1:10000), IP, etc.
Specificity	Reacts with NS4a polypeptide of Hepatitis C virus (a.a. 1055-1064) (RDKNQVEGEV).
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-004-018

