


# gp350/220 (RBD/(a.a.4-450))(EBV)



Product Name	gp350/220 (RBD/(a.a.4-450))(EBV)
Catalog No	IT-005-034p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged gp350/220 (RBD/(a.a.4-450))(EBV) (a.a.4-450) (See its amino acid sequence below; Genebank #: AHA36359). The receptor binding region of the gp350 protein.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



## Protein SEQ:

Allvcqytiqslihltgedpgffnveipefpyptcnvctadvnvtinfdvvgkqhldldfgqltphtkavyqprgafggse  
natnfllellgagelaltmrskklpinvttgeeqqvslesvdvyfqdvfgtmwchhaemqnpvylipetvpyikwdncnst  
nitavvraqgldvtlplsptsaqdsnfsvktemlgneidiecedgeisqvlpgdnkfnitcsgyeshvpsggiltstspvat  
pipgtgyaysrltprpvsrflgnnsilyvfysngpkasggdyciqsnivfsdeipasqdmptnttdityvgdnatysvpmv  
tsedanspnvtvtafwawpnttetdfkckwltlsgtspgcnisgafasnrtfditvsglgtapktliitrtatnattthkviska  
pestttsptlnttgfadpnttglps