


NS4 (HCV/Genotype-1a/Con)



Product Name	NS4 (HCV/Genotype-1a/Con)
Catalog No	IT-004-094Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein (1 mg/ml) in PBS with 8M urea.
Protein	GST-fused HCV NS4 Genotype-1a consensus immunodominant regions (aa 1659-1854).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Stable for 1 year at -80°C from the date of shipment. Avoiding repeated freeze/thaw. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

stwvlvggvlaalaayclstgcvvivrivlsgkpaipdrevlyrefdemeeesqhlpyieqgmmlaeqfkqkalgllqtas
rqaeviavqtnwqkleafwakhmwnfisiqylaglstlpgnpaiaslmaftaavtspltsqtlfnilggwvaaqlaap
gaatafvgaglagaaigsvlgkvlvdila