



gp120(MG505.W0M.ENV.H3)

Product Name	gp120(MG505.W0M.ENV.H3)(Clade A)
Catalog No	IT-001-102p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (MG505.W0M.ENV.H3)(Clade A) protein (a.a.31-508) (Genebank No.DQ208455).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E

E: gp120(MG505.W0M.ENV.H3) protein

Note: