



# gp120(MK184.W0M,ENV.E4)

---

Product Name	gp120(MK184.W0M,ENV.E4)(Clade A)
Catalog No	IT-001-111p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HIV-1 gp120 (MK184.W0M,ENV.E4)(Clade A) protein (a.a.35-510) (Genebank No.DQ208485).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(MK184.W0M,ENV.E4) protein

Note: