

gp120(00KE_KER2018) (HIV-1/Clade A)



Product Name	gp120(00KE_KER2018) (HIV-1/Clade A)
Catalog No	IT-001-149p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp120 (00KE_KER2018) (HIV-1/Clade A) protein (a.a. 31-504) (Genebank No. AAN05642)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

