

gp120(HIV-1/Clade 01_AE/90TH_CM244)



Product Name	gp120(HIV-1/Clade 01_AE/90TH_CM244)
Catalog No	IT-001-152p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp120(HIV-1/Clade 01_AE/90TH_CM244) protein (a.a.30-503)(Genebank #: AY713425)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(HIV-1/Clade 01_AE/90TH_CM244)