

Biotinylated gp140(HIV-1/Clade B'/C(CN54))



| | |
|----------------|--|
| Product Name | Biotinylated gp140(HIV-1/Clade B'/C(CN54)) |
| Catalog No | IT-001-161p-Biotin |
| Source | Recombinant protein purified from 293 cell culture |
| Formulation | Each vial contains 50 µg of biotinylated protein (0.5 mg/ml) in PBS+40% glycerol |
| Protein | Biotin labeled C-terminal 6xHis tagged gp140 (HIV-1/Clade B'/C(CN54)) protein (Genebank #: AX149771) (a.a. 29-676) |
| Applications | antibody ELISA, protein labeling, etc |
| Purity | >95% |
| Storage | Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |