

# HA( $\Delta$ TM) (B/Victoria/2/87)



---

Product Name	HA( $\Delta$ TM) (B/Victoria/2/87)
Catalog No	IT-003-B6p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM) (B/Victoria/2/87) protein (a.a. 16-548) (Genebank No. AAA43697)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

