



HA₁(A/swine/Guangxi/13/2006(H₁N₂))

Product Name	HA ₁ (A/swine/Guangxi/13/2006 (H ₁ N ₂))
Catalog No	IT-003-SW9p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA ₁ (A/swine/Guangxi/13/2006 (H ₁ N ₂))protein (a.a.18-342) (Genebank No. ABQ42444)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

