

NS5a (HCV/Genotype-1b)



Product Name	NS5a (HCV/Genotype-1b)
Catalog No	IT-004-079aEp
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	N-terminal GST, C-terminal 6xHis tagged HCV NS5a Genotype1b immunodominant regions, amino acids 2212-2321
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

