


gp140(aBL01)(HIV-1/Clade B)



Product Name	gp120(aBL01)(HIV-1/Clade B)
Catalog No	IT-001-171p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged gp140 (aBL01)(HIV-1/Clade B) protein(a.a. 31-667)(GenBank #: AAN39728)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140 (aBL01)(HIV-1/Clade B)