

GP1 [LASSA Virus]



Product Name	GP1 [LASSA Virus]
Catalog No	IT-008-002p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged receptor binding domain Of Lassa Virus GP1 (aa 16-172; GenBank# 4ZJF_A).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	Cohen-Dvashi, H., Cohen, N., Israeli, H. and Diskin, R. (2015). Molecular Mechanism for LAMP1 Recognition by Lassa Virus. J. Virol. 89, 7584-92.

