


# Human RANK L/TRANCE/TNFSF11



Product Name	Human RANK L/TRANCE/TNFSF11
Catalog No	IT-300-094Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl, and 8M urea.
Protein	6xHis tagged Human RANK L/TRANCE/TNFSF11 (a.a. 140–317) (GenBank Accession No. O14788).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human RANK L/TRANCE/TNFSF11