


# NA(A/Puerto Rico/8/1934(H1N1))



Product Name	NA(A/Puerto Rico/8/1934(H1N1))
Catalog No	IT-003-1012p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged NA(A/Puerto Rico/8/1934(H1N1)) (aa 36-454; GenBank# ABD77678).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥80%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: NA(A/Puerto Rico/8/1934(H1N1))

Note: