




# Human Podoplanin / PDPN

Product Name	Human Podoplanin / PDPN
Catalog No	IT-300-171p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of purified protein (1 mg/ml) in PBS.
Protein	Human Podoplanin / PDPN (aa 23-131; GenBank# Q86YL7) with 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human Podoplanin / PDPN

Note: