

# gD (HSV-2)



---

Product Name	gD (HSV-2)
Catalog No	IT-005-056Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M urea.
Protein	The E.Coli derived 39.7 kDa recombinant protein contains the HSV-2 gD immunodominant region 266-393 amino acids (GenBank# AEV91405), fused with 26 kDa GST-tag.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>92% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

