


# Human Interleukin-6 / IL-6 / BSF-2 / IFN-beta-2



Product Name	Human Interleukin-6 / IL-6 / BSF-2 / IFN-beta-2
Catalog No	IT-300-183Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl and 8M urea.
Protein	Human Interleukin-6 / IL-6 / BSF-2 / IFN-beta-2 (a.a. 30 – 212) (GenBank Accession No. P05231) with N-terminal 6xHis tag.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human Interleukin-6 / IL-6 / BSF-2 / IFN-beta-2