


## gp120(ZM215F)(HIV/Clade C)



Product Name	gp120(ZM215F)(HIV/Clade C)
Catalog No	IT-001-RC18p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	gp120(ZM215F)(HIV/Clade C) (aa 31-479; GenBank# DQ422948) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: gp120(ZM215F)(HIV/Clade C)

Note: