



Human NKG2A/CD94 Complex

Product Name	Human NKG2A/CD94 Complex
Catalog No	IT-300-190p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human NKG2A/CD94 Complex (aa 113-233, Genbank No. P26715; aa 57-179, Genbank No. Q13241) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

