



HA(Δ TM)(H5N1)(A/feline/Texas/24-029329-02/2024)

Product Name	HA(Δ TM)(H5N1)(A/feline/Texas/24-029329-02/2024)
Catalog No	IT-003-00542 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	The HA(Δ TM)(H5N1)(A/feline/Texas/24-029329-02/2024)(aa 17-532, GenBank # WZL56519) with 6xHis tag at C-terminus. Majority in trimer.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	\geq 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H5N1)
(A/feline/Texas/24-029329-02/2024)

Note: