

# Anti- Matrix protein 1 [Dhori virus]



Product Name	Anti- Matrix protein 1 [Dhori virus]
Catalog No	IT-025-006M2
Source	Mouse monoclonal antibody, clone 2G8F7
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Recombinant Matrix protein 1 [Dhori virus] (IT-025-006Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Matrix protein 1 [Dhori virus] as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

